


Cadenza CD-C18

Technical

Reproducibility of batch-to-batch



This data shows Cadenza CD-C18's high level of repetition between packing gel production batches.

This experiment separated a drug-related chemical compound using three columns packed with ODS packings from different batches. The eluent is neutral pH and LC-MS compatible with volatile and low ionic strength.

Even in this quick analysis of multiple ingredients, Cadenza CD-C18's elution peaks show no change in elution peak behavior. This proves high-reproducibility even with multiple ingredients including base compounds.

Cadenza CD-C18's material processing, ODS substitution method, and end-capping technology guarantees a high level of repetition with its superior total design, manufacture and quality control techniques.