


Unison UK-Phenyl  
Unison UK-C8  
Cadenza CW-C18  
Cadenza CD-C18  
Unison UK-C18

150 x 4.6 mm

Application

## Cold medicine compounds



**Unison UK-Phenyl**, 150 x 4.6 mm  
100mM ammonium formate / acetonitrile = 80 / 20  
1 mL/min (9MPa), 37 deg.C, 235 nm, 2 uL (0.4-1.8 ug)

**Unison UK-C8**, 150 x 4.6 mm  
100mM ammonium formate / methanol = 60 / 40  
1 mL/min (16MPa), 37 deg.C, 235 nm, 2 uL (0.4-1.8 ug)

**Cadenza CW-C18**, 150 x 4.6 mm  
100mM ammonium formate / methanol = 60 / 40  
1 mL/min (17MPa), 37 deg.C, 235 nm, 2 uL (0.4-1.8 ug)

**Unison UK-C18**, 150 x 4.6 mm  
100mM ammonium formate / methanol = 50 / 50  
1 mL/min (16MPa), 37 deg.C, 235 nm, 2 uL (0.4-1.8 ug)

**Cadenza CD-C18**, 150 x 4.6 mm  
100mM ammonium formate / methanol = 50 / 50  
1 mL/min (17MPa), 37 deg.C, 235 nm, 2 uL (0.4-1.8 ug)