

**Intrada Amino Acid**

50 x 3 mm

Application

**55 Amino acids (standard)**

55種アミノ酸(標準試料)

**Intrada Amino Acid, 50 x 3 mm**

**A: acetonitrile /tetrahydrofuran /25mM ammonium formate /formic acid = 9 / 75 / 16 / 0.3**

**B: acetonitrile / 100mM ammonium formate = 20 / 80**

**0 %B (0-2.5min)**

**0-17 %B (2.5-6.5 min)**

**100 %B (6.5-10 min)**

**0.6 mL/min (6MPa), 35 deg.C, 5 uL (1000nmol/mL)**

**Shimadzu LCMS-2020, ESI (SIM, positive)**

